


Caveolin-1 Polyclonal Antibody (Rabbit) For Wb, If-P, Ihc-P, Elisa - Q03135

Item Number: CM0023732



Introduction

Rabbit polyclonal antibody targeting human Caveolin-1 (CAV1), validated for Western blot (WB), immunofluorescence (IF-P), immunohistochemistry (IHC-P), and ELISA. Cross-reacts with mouse and rat. Suitable for caveolae research, cell signaling, and membrane trafficking studies.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-1 Rabbit pAb
Remarks/Alias	CGL3; PPH3; BSCL3; LCCNS; VIP21; MSTP085; Caveolin-1
Species	Human
Gene ID (Human)	857
Gene ID	857
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human Caveolin-1 (NP_001744.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q03135
Protein Weight	20kDa
Shipping	Ice bag