

Caveolin-3 Rabbit Pab - P56539

Item Number: CM0016008

Introduction

Caveolin-3 Rabbit pAb for WB and ELISA, validated in human, mouse, rat. Targets caveolin-3, a muscle-specific scaffolding protein for membrane repair and signaling.

[Learn More](#)

Parameter	Value
Product Name	Caveolin-3 Rabbit pAb
Remarks/Alias	LQT9; MPDT; RMD2; VIP21; LGMD1C; VIP-21; Caveolin-3
Species	Human
Gene ID	859
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Caveolin-3 (NP_001225.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P56539
Protein Weight	17kDa
Shipping	Ice bag