

Cbr3 (Sdr21C2) Rabbit Polyclonal Antibody For Wb, Elisa - 075828

Item Number: CM0024769



Introduction

High-quality CBR3 rabbit polyclonal antibody validated for WB and ELISA. Recognizes human, mouse, rat CBR3. Targets carbonyl reductase 3, involved in NADPH-dependent reduction of quinones and anticancer drugs. Ideal for IVD research and development.

[Learn More](#)

Parameter	Value
Product Name	CBR3 Rabbit pAb
Remarks/Alias	hCBR3; SDR21C2; HEL-S-25; CBR3
Species	Human
Gene ID (Human)	874
Gene ID	874
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 118-277 of human CBR3 (NP_001227.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	075828
Protein Weight	31kDa
Shipping	Ice bag