


Cbx2 Rabbit Monoclonal Antibody For Wb & Elisa - Q14781

Item Number: CM0011690



Introduction

Rabbit monoclonal antibody to human CBX2, a Polycomb group protein. Reacts with human, mouse, rat. Validated for WB and ELISA. CBX2 is critical for chromatin regulation, gene silencing, and sexual development. Ideal for epigenetic research.

[Learn More](#)

Parameter	Value
Product Name	CBX2 Rabbit mAb
Remarks/Alias	M33; CDCA6; SRXY5; CBX2
Species	Human
Gene ID	84733
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 433-532 of human CBX2 (Q14781)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14781
Protein Weight	56kDa
Shipping	Ice bag