

Ccar2/Dbc1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q8N163

Item Number: CM0019490



Introduction

Rabbit polyclonal antibody against CCAR2/DBC1 (p30 DBC) for WB, IHC-P, ELISA. Detects human, mouse, rat protein. Investigate roles in apoptosis, transcriptional regulation, and circadian rhythms.

[Learn More](#)

Parameter	Value
Product Name	CCAR2/DBC1 Rabbit pAb
Remarks/Alias	CCAR2; DBC-1; DBC1; KIAA1967; NET35; p30 DBC; p30DBC; cell cycle and apoptosis regulator 2; CCAR2/DBC1
Species	Human
GeneID (Human)	57805
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CCAR2/DBC1 (NP_066997.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N163
Protein Weight	102kDa/103kDa
Shipping	Ice bag