


Ccl21 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - 000585

Item Number: CM0021734



Introduction

Rabbit polyclonal antibody targeting human CCL21 (6Ckine/SLC), validated for WB, IHC-P, and ELISA. Ideal for lymphocyte homing and chemokine research. Recognizes a band at ~15 kDa. Reacts with human CCL21.

[Learn More](#)

Parameter	Value
Product Name	CCL21 Rabbit pAb
Remarks/Alias	ECL, SLC, CKb9, TCA4, 6Ckine, SCYA21, CCL21
Species	Human
GeneID	6366
Immunogen	Recombinant fusion protein containing amino acids 24-134 of human CCL21 (NP_002980.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O00585
Protein Weight	15kDa
Shipping	Ice bag