

Ccl25 Rabbit Pab For Wb, Elisa - O15444

Item Number: CM0019213

Introduction

CCL25 Rabbit pAb is a polyclonal antibody specific to human CCL25, validated for Western blot and ELISA. Ideal for studying thymus-expressed chemokine functions in T-cell development and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	CCL25 Rabbit pAb
Remarks/Alias	TECK; Ckb15; SCYA25; TECKvar; Ck beta-15; CCL25
Species	Human
GeneID (Human)	6370
GeneID	6370
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-150 of human CCL25 (NP_005615.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15444
Protein Weight	17kDa
Shipping	Ice bag