

Ccl7 Rabbit Pab For Wb, Elisa - P80098

Item Number: CM0022573



Introduction

Rabbit polyclonal CCL7 antibody validated for Western blotting and ELISA. Reacts with human and rat CCL7, a chemokine involved in immune cell recruitment and inflammation. Ideal for immunology and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CCL7 Rabbit pAb
Remarks/Alias	FIC; MARC; MCP3; NC28; MCP-3; SCYA6; SCYA7; CCL7
Species	Human
GeneID (human)	6354
GeneID	6354
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-99 of human CCL7 (NP_006264.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P80098
Protein Weight	11kDa
Shipping	Ice bag