

Ccr7 Rabbit Pab - P47774

Item Number: CM0030757

Introduction

CCR7 Rabbit Polyclonal Antibody for WB and ELISA. Reacts with mouse and rat CCR7. Suitable for investigating chemokine receptor-mediated immune cell trafficking and inflammatory responses.

[Learn More](#)

Parameter	Value
Product Name	CCR7 Rabbit pAb
Remarks/Alias	EBI1; CCR-7; CD197; Ebi1h; Cdw197; Cmkbr7; CC-CKR-7
Species	Mouse
GeneID (Human)	1236
GeneID	12775
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 211-310 of mouse CCR7 (NP_031745.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P47774
Protein Weight	43kDa
Shipping	Ice bag