

Cd11B (Integrin Alpha-M) Rabbit Polyclonal Antibody For Wb And Elisa - P11215

Item Number: CM0019092



Introduction

Rabbit polyclonal antibody targeting human CD11b (Integrin alpha-M), validated for Western blot and ELISA applications with cross-reactivity to mouse and rat. Suitable for immune cell adhesion, complement receptor, and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD11b Rabbit pAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6; CD11b/ITGAM
Species	Human
Gene ID (Human)	3684
Gene ID	3684
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD11b (NP_000623.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag