

Cd11B Rabbit Polymab® Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11215 / P05555

Item Number: CM0009373



Introduction

Rabbit monoclonal antibody targeting human CD11b (Integrin α -M, MAC-1) validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for studies on innate immunity, cell adhesion, and leukocyte function.

[Learn More](#)

Parameter	Value
Product Name	CD11b Rabbit PolymAb®
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6
Species	Human
Gene ID (Human)	3684
Gene ID	3684/16409
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human CD11b (NP_000623.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11215/P05555
Protein Weight	127kDa
Shipping	Ice bag