

Cd11B/Itgam Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa - P11215

Item Number: CM0014983



Introduction

Rabbit polyclonal antibody targeting integrin alpha-M (CD11b/ITGAM). Validated for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with human, mouse, rat. Suitable for studying cell adhesion, phagocytosis, and innate immunity.

[Learn More](#)

Parameter	Value
Product Name	CD11b/ITGAM Rabbit pAb
Remarks/Alias	CR3A; MO1A; CD11B; MAC-1; MAC1A; SLEB6; CD11b/ITGAM
Species	Human
GeneID	3684 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 145-340 of CD11b/ITGAM (NP_000623.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11215
Protein Weight	127kDa
Shipping	Ice bag