

Cd138 Rabbit Mab - P18827

Item Number: CM0008227



Introduction

CD138 Rabbit monoclonal antibody for IHC-P, FC, and ELISA. Recognizes human and cynomolgus monkey CD138/Syndecan-1. Ideal for plasma cell research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	CD138 Rabbit mAb
Remarks/Alias	SDC; CD138; SYND1; syndecan
Species	Human
Gene ID (Human)	6382
Gene ID	6382
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-251 of human CD138 (NP_002988.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, FC, ELISA
Cross Reactivity	Human, Cynomolgus monkey
SWISS	P18827
Protein Weight	32kDa
Shipping	Ice bag