

Cd14 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08571

Item Number: CM0028420



Introduction

CD14 rabbit polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects CD14 in human, mouse, and rat samples. Ideal for innate immunity and LPS signaling research.

[Learn More](#)

Parameter	Value
Product Name	CD14 Rabbit pAb
Remarks/Alias	CD14
Species	Human
GeneID (Human)	929
GeneID	929
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-330 of human CD14 (NP_000582.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08571
Protein Weight	40kDa
Shipping	Ice bag