


Cd146/Mcam Rabbit Polyclonal Antibody For Wb & Elisa - P43121

Item Number: CM0021970



Introduction

High-quality rabbit polyclonal antibody targeting human CD146/MCAM, validated for Western blotting and ELISA. Detects human and rat MCAM, ideal for cell adhesion, cancer metastasis, and vascular biology research.

[Learn More](#)

Parameter	Value
Product Name	CD146/MCAM Rabbit pAb
Remarks/Alias	CD146, MUC18, HEMCAM, METCAM, MelCAM
Species	Human
GeneID	4162
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 447-646 of human CD146/MCAM (NP_006491.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P43121
Protein Weight	72kDa
Shipping	Ice bag