

Cd160 Rabbit Pab For Wb And Elisa - O95971

Item Number: CM0028577



Introduction

Rabbit polyclonal antibody targeting human CD160, validated for WB and ELISA, with cross-reactivity to mouse and rat. Ideal for immune cell receptor studies and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD160 Rabbit pAb
Remarks/Alias	NK1, BY55, NK28, CD160
Species	Human
Gene ID (Human)	11126
Gene ID	11126
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of human CD160 (NP_008984.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O95971
Protein Weight	20kDa
Shipping	Ice bag