

## Cd172A/Sirpα Rabbit Mab For Wb, If/Icc, Elisa - P78324

Item Number: CM0005496



### Introduction

CD172a/SIRPα Rabbit mAb (P78324) validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. High specificity monoclonal antibody for signal-regulatory protein alpha research.

[Learn More](#)

| Parameter        | Value                                                                                                                 |
|------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name     | CD172a/SIRPα Rabbit mAb                                                                                               |
| Remarks/Alias    | BIT; MFR; P84; SIRP; MYD-1; SHPS1; CD172A; PTPNS1; CD172a/SIRPα                                                       |
| Species          | Human                                                                                                                 |
| Gene ID          | 140885                                                                                                                |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 395-495 of human CD172a/SIRPα (P78324). |
| Source           | Rabbit                                                                                                                |
| Category         | Monoclonal Antibodies                                                                                                 |
| Application      | WB, IF/ICC, ELISA                                                                                                     |
| Cross Reactivity | Human, Mouse, Rat                                                                                                     |
| SWISS            | P78324                                                                                                                |
| Protein Weight   | 55kDa                                                                                                                 |
| Shipping         | Ice bag                                                                                                               |