

Cd18 Rabbit Pab - P05107

Item Number: CM0024283



Introduction

CD18 Rabbit Polyclonal Antibody detects human Integrin beta-2 (ITGB2, LFA-1, MAC-1) and cross-reacts with mouse and rat. Validated for WB, IHC-P, IF/ICC, ELISA, ideal for leukocyte adhesion and immune response research.

[Learn More](#)

Parameter	Value
Product Name	CD18 Rabbit pAb
Remarks/Alias	LAD; CD18; MF17; MFI7; LCAMB; LFA-1; MAC-1
Species	Human
GeneID (Human)	3689
GeneID	3689
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-125 of CD18 (NP_000202.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05107
Protein Weight	85kDa
Shipping	Ice bag