

Cd19 Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P15391

Item Number: CM0003393



Introduction

CD19 rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with human and cynomolgus monkey CD19, a B-cell coreceptor involved in BCR signaling. Ideal for B-cell research.

[Learn More](#)

Parameter	Value
Product Name	CD19 Rabbit mAb
Remarks/Alias	B4; CVID3; CD19
Species	Human
Gene ID	930
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 457-556 of human CD19 (NP_001761.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Cynomolgus monkey
SWISS	P15391
Protein Weight	61kDa
Shipping	Ice bag