


Cd1B Rabbit Pab For Wb, Elisa - P29016

Item Number: CM0030016



Introduction

Rabbit polyclonal antibody targeting human CD1B, a T-cell surface glycoprotein that presents lipid antigens to NKT cells. Suitable for Western blot and ELISA. Validated with recombinant immunogen.

[Learn More](#)

Parameter	Value
Product Name	CD1B Rabbit pAb
Remarks/Alias	R1; CD1; CD1A; CD1B
Species	Human
Gene ID (Human)	910
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-303 of human CD1B (NP_001755.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P29016
Protein Weight	37kDa
Shipping	Ice bag