


Cd2 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P06729

Item Number: CM0029236



Introduction

High-quality rabbit polyclonal antibody against human CD2 (T11, LFA-2). Validated for WB, IHC-P, and ELISA. Ideal for studying T-cell adhesion, lymphocyte activation, and cell signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	CD2 Rabbit pAb
Remarks/Alias	T11; SRBC; LFA-2; CD2
Species	Human
GeneID (Human)	914
GeneID	914
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 235-351 of human CD2 (NP_001758.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P06729
Protein Weight	39kDa
Shipping	Ice bag