

Cd267/Taci Rabbit Monoclonal Antibody For Flow Cytometry - Q9Et35

Item Number: CM0006252

Introduction

Recombinant rabbit monoclonal antibody targeting mouse CD267/TACI (TNFRSF13B), validated for flow cytometry. Receptor for APRIL and BAFF, key in B-cell function and humoral immunity. For immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD267/TACI Rabbit mAb
Remarks/Alias	Taci; 1200009E08Rik; CD267/TACI
Species	Mouse
GeneID (Human)	23495
GeneID	57916
Immunogen	Recombinant Protein Recombinant fusion protein containing a sequence corresponding to amino acids 5-129 of mouse CD267/TACI (NP_067324.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Mouse
SWISS	Q9ET35
Protein Weight	27kDa
Shipping	Ice bag