

Cd275 Rabbit Mab For Wb, Elisa - O75144

Item Number: CM0005238



Introduction

High-quality rabbit monoclonal antibody targeting human CD275 (B7H2/ICOSL). Validated for WB and ELISA applications, ideal for T-cell costimulation and immune response studies.

[Learn More](#)

Parameter	Value
Product Name	CD275 Rabbit mAb
Remarks/Alias	B7h; B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L
Species	Human
GeneID (Human)	23308
GeneID	23308
Immunogen	Recombinant protein of human CD275.
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75144
Protein Weight	33kDa
Shipping	Ice bag