

Cd31/Pecam1 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P16284

Item Number: CM0019027



Introduction

Rabbit polyclonal antibody against CD31/PECAM1, validated for western blot, immunofluorescence, and ELISA. Reacts with human, mouse, and rat. Recognizes an 83kDa protein critical for leukocyte transendothelial migration, cell adhesion, and phagocytosis. Suitable for cardiovascular and inflammation research.

[Learn More](#)

Parameter	Value
Product Name	CD31/PECAM1 Rabbit pAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM; CD31/PECAM1
Species	Human
GeneID (Human)	5175
GeneID	5175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 630-738 of CD31/PECAM1 (NP_000433.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P16284
Protein Weight	83kDa
Shipping	Ice bag