

Cd34 Rabbit Pab For Wb - P28906

Item Number: CM0014361

Introduction

CD34 Rabbit Polyclonal Antibody, validated for WB, IF-P, IHC-P, and ELISA in human and rat. Suitable for hematopoiesis and stem cell adhesion studies. 41kDa.

[Learn More](#)

Parameter	Value
Product Name	CD34 Rabbit pAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID (Human)	947
GeneID	947
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 32-290 of human CD34 (NP_001020280.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Rat
SWISS	P28906
Protein Weight	41kDa
Shipping	Ice bag