

Cd4 Rabbit Polyclonal Antibody For Wb, If-P, Elisa - P01730

Item Number: CM0015716



Introduction

CD4 rabbit polyclonal antibody validated for WB, IF-P, ELISA in human samples. Recognizes T-cell surface glycoprotein CD4, a critical immune coreceptor and HIV receptor. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	CD4 Rabbit pAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut; CD4
Species	Human
GeneID (human)	920
GeneID	920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-207 of CD4 (NP_000607.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P01730
Protein Weight	51 kDa
Shipping	Ice bag