

Cd42B Rabbit Mab For Wb, Ihc-P, Elisa - P07359

Item Number: CM0006872



Introduction

CD42b Rabbit mAb (SWISS: P07359) validated for WB, IHC-P, and ELISA. Targets human platelet glycoprotein Ib alpha chain involved in hemostasis. Ideal for platelet biology and bleeding disorder research.

[Learn More](#)

Parameter	Value
Product Name	CD42b Rabbit mAb
Remarks/Alias	BSS; GP1B; VWDP; CD42B; GPIbA; BDPLT1; BDPLT3; DBPLT3; GPIbalpha; CD42b-alpha; CD42b
Species	Human
GeneID (Human)	2811
GeneID	2811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-531 of human CD42b (NP_000164.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P07359
Protein Weight	72kDa
Shipping	Ice bag