

Cd43 Rabbit Mab For Wb, If/Icc, Elisa - P16150

Item Number: CM0006971

Introduction

Monoclonal rabbit anti-CD43 antibody validated for western blot, IF/ICC, and ELISA on human samples. Targets leukosialin, a key regulator of T-cell functions and leukocyte recruitment.

[Learn More](#)

Parameter	Value
Product Name	CD43 Rabbit mAb
Remarks/Alias	SPN; CD43; GALGP; GPL115; LSN; leukosialin
Species	Human
Gene ID	6693
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-253 of human CD43 (NP_003114.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P16150
Protein Weight	40kDa
Shipping	Ice bag