


Cd44 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P16070

Item Number: CM0014261



Introduction

Rabbit polyclonal antibody against human CD44, validated for Western blot, immunohistochemistry (paraffin), and ELISA. Cross-reacts with mouse. Suitable for cell adhesion, migration, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	CD44 Rabbit pAb
Remarks/Alias	IN, LHR, MC56, MDU2, MDU3, MIC4, Pgp1, CDW44, CSPG8, H-CAM, HCELL, ECM-III, HUTCH-1, HUTCH-I, ECMR-III, Hermes-1, CD44
Species	Human
GeneID (Human)	960
GeneID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of CD44 (NP_000601.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag