

Cd53 Polyclonal Antibody (Rabbit) For Wb And Elisa - Q61451

Item Number: CM0023622

Introduction

Rabbit polyclonal antibody against mouse CD53 (Ox-44, Tspan25) for WB and ELISA. Validated in mouse samples; ideal for studying tetraspanin functions in B-cell development, lymphocyte homing, and inflammatory signaling.

[Learn More](#)

Parameter	Value
Product Name	CD53 Rabbit pAb
Remarks/Alias	Ox-44; Tspan25; CD53
Species	Mouse
GeneID	12508
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 107-181 of mouse CD53 (NP_031677.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q61451
Protein Weight	24kDa
Shipping	Ice bag