

Cd7 Rabbit Pab - P09564

Item Number: CM0026805

Introduction

CD7 Rabbit Polyclonal Antibody validated for WB, FC, ELISA. Cross-reacts with human and rat T-cell antigen CD7 (GP40/TP41). Ideal for immunology research on T-cell and NK cell activation.

[Learn More](#)

Parameter	Value
Product Name	CD7 Rabbit pAb
Remarks/Alias	GP40, TP41, Tp40, LEU-9, CD7
Species	Human
GeneID (Human)	924
GeneID	924
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, FC, ELISA
Cross Reactivity	Human, Rat
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag