

Cd79A Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Fc, Elisa - P11912

Item Number: CM0007249



Introduction

High-quality CD79a Rabbit Monoclonal Antibody validated for Human samples in Western blot, IHC-P, IF/ICC, flow cytometry, and ELISA.

Detects the 25kDa B-cell receptor alpha chain (SWISS P11912), crucial for BCR signaling and immune research.

[Learn More](#)

Parameter	Value
Product Name	CD79a Rabbit mAb
Remarks/Alias	IGA; MB1; MB-1; IGA1pha; CD79a
Species	Human
GeneID (Human)	973
GeneID	973
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 165-225 of human CD79a (NP_001774.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, FC, ELISA
Cross Reactivity	Human
SWISS	P11912
Protein Weight	25kDa
Shipping	Ice bag