

Cd79B Rabbit Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P40259

Item Number: CM0004829



Introduction

CD79B Rabbit Monoclonal Antibody validated for WB, IF-P, IHC-P, ELISA. Reacts with human, mouse, rat. Suitable for B-cell receptor signaling research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	CD79B Rabbit mAb
Remarks/Alias	B29; IGB; AGM6; Igbeta; CD79B
Species	Human
GeneID (Human)	974
GeneID	974
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P40259
Protein Weight	26kDa
Shipping	Ice bag