

Cdc25C Rabbit Pab For If/Icc, Elisa - P30307

Item Number: CM0017004



Introduction

CDC25C Rabbit pAb is a polyclonal antibody validated for IF/ICC and ELISA, targeting the M-phase inducer phosphatase 3. Suitable for human samples in cell cycle and oncology research.

[Learn More](#)

Parameter	Value
Product Name	CDC25C Rabbit pAb
Remarks/Alias	CDC25; PPP1R60; CDC25C
Species	Human
GeneID (Human)	995
GeneID	995
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P30307
Protein Weight	53kDa
Shipping	Ice bag