

Cdc5L Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q99459

Item Number: CM0011029



Introduction

CDC5L Rabbit mAb, validated for WB, IHC-P, ELISA. Recognizes human, mouse CDC5L (92 kDa). Ideal for splicing, cell cycle, DNA damage studies. UniProt Q99459.

[Learn More](#)

Parameter	Value
Product Name	CDC5L Rabbit mAb
Remarks/Alias	CDC5; CEF1; PCDC5RP; CDC5-LIKE; dj319D22.1; CDC5L
Species	Human
GeneID (Human)	988
GeneID	988
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99459
Protein Weight	92kDa
Shipping	Ice bag