

Cdh9 Rabbit Polyclonal Antibody For Wb And Elisa - Q9ULb4

Item Number: CM0013028

Introduction

High-quality rabbit polyclonal antibody against human CDH9 (Cadherin-9), validated for WB and ELISA. Ideal for cell adhesion research. Cross-reacts with human. UniProt: Q9ULB4.

[Learn More](#)

Parameter	Value
Product Name	CDH9 Rabbit pAb
Remarks/Alias	CDH9
Species	Human
GeneID (Human)	1007
GeneID	1007
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-420 of human CDH9 (NP_057363.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9ULB4
Protein Weight	89kDa
Shipping	Ice bag