


Cdk1 Rabbit Monoclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P06493

Item Number: CM0005064



Introduction

CDK1 Rabbit monoclonal antibody validated for WB, IHC-P, IF/ICC, and ELISA. Detects CDK1 in human, mouse, and rat samples. CDK1 is a critical serine/threonine kinase driving G2/M transition and mitosis. Ideal for oncology, cell cycle, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	CDK1 Rabbit mAb
Remarks/Alias	CDC2; CDC28A; P34CDC2; CDK1
Species	Human
GenelD	983
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P06493
Protein Weight	34kDa
Shipping	Ice bag