

Cdk4 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P11802

Item Number: CM0006590



Introduction

High-quality rabbit monoclonal antibody targeting human CDK4 (PSK-J3). Validated for WB, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Detects a 34kDa protein. Ideal for cell cycle, cancer, and RB pathway studies.

[Learn More](#)

Parameter	Value
Product Name	CDK4 Rabbit mAb
Remarks/Alias	CMM3; PSK-J3; CDK4
Species	Human
GeneID (Human)	1019
GeneID	1019
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 200-303 of human CDK4 (NP_000066.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11802
Protein Weight	34kDa
Shipping	Ice bag