

Cdk4 Rabbit Polyclonal Antibody [Ko Validated] For Wb, Ihc-P, If/Icc, Ip, Elisa - P11802

Item Number: CM0020958



Introduction

Rabbit polyclonal CDK4 antibody, KO validated, for WB, IHC-P, IF/ICC, IP, ELISA. Detects human, mouse, rat CDK4. Synthetic peptide immunogen (aa 204-303). Ideal for cell cycle and cancer research.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CDK4 Rabbit pAb
Remarks/Alias	CMM3; PSK-J3; K4
Species	Human
GeneID (Human)	1019
GeneID	1019
Immunogen	Synthetic peptide corresponding to amino acids 204-303 of human CDK4 (NP_000066.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11802
Protein Weight	34kDa
Shipping	Ice bag