

Cdk9 Rabbit Polyclonal Antibody For Wb And Elisa - P50750

Item Number: CM0013116



Introduction

Rabbit polyclonal antibody against CDK9, validated for Western blot (WB) and ELISA. Reacts with human, mouse, and rat. Suitable for studies of transcription regulation, cell cycle, and chromatin modification.

[Learn More](#)

Parameter	Value
Product Name	CDK9 Rabbit pAb
Remarks/Alias	TAK; C-2k; CTK1; CDC2L4; PITALRE; CDK9
Species	Human
GeneID (Human)	1025
GeneID	1025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 153-372 of human CDK9 (NP_001252.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P50750
Protein Weight	43kDa
Shipping	Ice bag