


Cdkn2A/P16Ink4A Rabbit Pab For Wb, Elisa - P42771 / Q8N726

Item Number: CM0017794



Introduction

Rabbit polyclonal antibody targeting human CDKN2A/p16INK4a, validated for WB and ELISA. Suitable for cell cycle and tumor suppressor studies. Recognizes multiple isoforms (8-17 kDa). Cross-reacts with human.

[Learn More](#)

Parameter	Value
Product Name	CDKN2A/p16INK4a Rabbit pAb
Remarks/Alias	ARF; MLM; P14; P16; P19; CMM2; INK4; MTS1; TP16; CDK4I; CDKN2; INK4A; MTS-1; P14ARF; P19ARF; P16INK4; P16INK4A; P16-INK4A; CDKN2A/p16INK4a
Species	Human
GeneID (Human)	1029
GeneID	1029
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human CDKN2A/p16INK4a (NP_000068.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P42771/Q8N726
Protein Weight	8kDa/11kDa/12kDa/13kDa/16kDa/17kDa
Shipping	Ice bag