

Cdv3 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uky7

Item Number: CM0020070

Introduction

CDV3 Rabbit Polyclonal Antibody suitable for WB, IF/ICC, and ELISA. Detects human, mouse, and rat CDV3 protein. Raised against a synthetic peptide within amino acids 150-250 of human CDV3. For research use in breast cancer and cell biology.

[Learn More](#)

Parameter	Value
Product Name	CDV3 Rabbit pAb
Remarks/Alias	H41; CDV3
Species	Human
Gene ID (Human)	55573
Gene ID	55573
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 150-250 of human CDV3 (NP_060018.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9UKY7
Protein Weight	27kDa
Shipping	Ice bag