

Cdx1 Rabbit Polyclonal Antibody For Wb, Elisa - P47902

Item Number: CM0026222



Introduction

CDX1 rabbit polyclonal antibody validated for Western blot (WB) and ELISA with human and mouse reactivity. Ideal for investigating transcriptional regulation in intestinal epithelium and colorectal cancer models.

[Learn More](#)

Parameter	Value
Product Name	CDX1 Rabbit pAb
Remarks/Alias	CDX1
Species	Human
GeneID (Human)	1044
GeneID	1044
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 136-265 of human CDX1 (NP_001795.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P47902
Protein Weight	28kDa
Shipping	Ice bag