

Cdx2 Rabbit Polyclonal Antibody For Ihc-P And Elisa - Q99626

Item Number: CM0019419



Introduction

CDX2 rabbit polyclonal antibody for IHC-P and ELISA. Reacts with human, mouse, rat CDX2 transcription factor. Ideal for intestinal development and differentiation studies.

[Learn More](#)

Parameter	Value
Product Name	CDX2 Rabbit pAb
Remarks/Alias	CDX3; CDX-3; CDX2/AS; CDX2
Species	Human
Gene ID (Human)	1045
Gene ID	1045
Immunogen	Recombinant Protein corresponding to a sequence within amino acids 1-180 of human CDX2 (NP_001256.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q99626
Protein Weight	34kDa
Shipping	Ice bag