

Ceacam7 Rabbit Polyclonal Antibody For Wb And Elisa - Q14002

Item Number: CM0025183



Introduction

Rabbit polyclonal antibody targeting CEACAM7 (CGM2) for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat. Ideal for studying cell adhesion in colon epithelium and colorectal cancer.

[Learn More](#)

Parameter	Value
Product Name	CEACAM7 Rabbit pAb
Remarks/Alias	CGM2; CEACAM7
Species	Human
GeneID (Human)	1087
GeneID	1087
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 36-242 of human CEACAM7 (NP_008821.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14002
Protein Weight	29kDa
Shipping	Ice bag