

Cend1 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - Q8N111

Item Number: CM0012112

Introduction

CEND1 Rabbit Monoclonal Antibody validated for WB, IF-P, ELISA. Reacts with human, mouse, rat. Suitable for neuronal differentiation research. UniProt Q8N111.

[Learn More](#)

Parameter	Value
Product Name	CEND1 Rabbit mAb
Remarks/Alias	BM88; CEND1
Species	Human
GeneID (Human)	51286
GeneID	51286
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 70-149 of human CEND1 (NP_057648.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N111
Protein Weight	15kDa
Shipping	Ice bag