

Cenpe Rabbit Polyclonal Antibody For Wb, Elisa - Q02224

Item Number: CM0022388



Introduction

CENPE Rabbit Polyclonal Antibody validated for WB and ELISA, reactive across human, mouse, and rat. Ideal for investigating chromosome congression and mitotic checkpoint signaling.

[Learn More](#)

Parameter	Value
Product Name	CENPE Rabbit pAb
Remarks/Alias	KIF10; CENP-E; MCPH13; PPP1R61; CENPE
Species	Human
GeneID (Human)	1062
GeneID	1062
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 2542-2701 of human CENPE (NP_001804.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q02224
Protein Weight	316kDa
Shipping	Ice bag