

Cep290 Rabbit Polyclonal Antibody For Wb, Elisa - O15078

Item Number: CM0024689



Introduction

High-quality rabbit polyclonal antibody targeting human CEP290 (290 kDa), a key centrosomal protein in cilia formation. Validated for Western blot and ELISA, supporting research into ciliopathies and centrosome biology.

[Learn More](#)

Parameter	Value
Product Name	CEP290 Rabbit pAb
Remarks/Alias	CT87; MKS4; POC3; rd16; BBS14; JBTS5; LCA10; NPHP6; SLSN6; 3H11Ag; CEP290
Species	Human
Gene ID (Human)	80184
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 950-1050 of human CEP290 (NP_079390.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O15078
Protein Weight	290kDa
Shipping	Ice bag