

Ces3 Rabbit Pab - Q6Uww8

Item Number: CM0029225

Introduction

High-quality CES3 rabbit polyclonal antibody validated for WB and ELISA applications. Detects endogenous CES3 in human, mouse, and rat samples. Ideal for studying xenobiotic detoxification and prodrug activation in liver, colon, and small intestine.

[Learn More](#)

Parameter	Value
Product Name	CES3 Rabbit pAb
Remarks/Alias	ES31; CES3
Species	Human
GeneID (Human)	23491
GeneID	23491
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 269-568 of human CES3 (NP_001172106.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6UWW8
Protein Weight	62kDa
Shipping	Ice bag