

Cga Rabbit Polyclonal Antibody For Wb And Elisa - P01215

Item Number: CM0017359



Introduction

Rabbit polyclonal antibody against human CGA (glycoprotein hormones alpha chain; ~13 kDa). Validated for Western blot and ELISA. Detects the shared alpha subunit of TSH, LH, FSH, and CG for endocrine and reproductive research.

[Learn More](#)

Parameter	Value
Product Name	CGA Rabbit pAb
Remarks/Alias	HCG; LHA; FSHA; GPA1; GPHa; TSHA; GPHA1; CG-ALPHA; CGA
Species	Human
GeneID (human)	1081
GeneID	1081
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P01215
Protein Weight	13kDa
Shipping	Ice bag